



# utopia

P R O T E C T S

## PRODUCT INFORMATION



Product - V30000  
KN95 Respirator Mask

Colour  
White

Material  
Outer: Non-woven material  
Middle: Two Layers of Melt-Blown material  
Inner: Non-woven material  
Adjustable nose clip

Quality Standard  
GB2626 2006 (EN149-2001 Equivalent)

Size  
One Size

Pack Size  
Box of 50, Individually wrapped

Instructions:  
With sanitised hands unfold the mask by pulling the mask apart. Put the mask over your nose and mouth. Adjust the straps to fit comfortably around the ears. Mould the metal bridge around your nose without pinching to make it fit snugly. Perform a seal check by covering the mask with both hands without pinching and inhale and exhale. If you feel air at the edge of the mask readjust it. Correct fit requires smooth contact with skin.

Discard after use  
Product life is two years from date of manufacture



*Serve and Protect*